



SILPEROX

miraculous disinfectant: SILPEROX
harmonic dance of
hydrogen peroxide & colloidal silver
(nano metallic silver)

What is SILPEROX




The most recent technological invention in disinfection area

based on patented - synergistic blend of

hydrogen peroxide and colloidal (non-ionic) silver,

thanks to the high nano-technology



What is SILPEROX



Thus,

germicidal potency,
cleaning performance,
and stability
characteristics dramatically increase.


Ensures the best disinfecting and cleaning performance




What is SILPEROX



effective against bacteria, viruses, fungi, yeasts, spores and even mycobacteria



kills a broad spectrum of pathogens in very short time



disinfects before drying on a surface

one-stop biocidal disinfectant product

What is SILPEROX




odorless,

chlorine free,



alcohol free,


and colorless



What is SILPEROX



does not irritate the skin and eyes (when it is used at recommended dosages)



Effective at high temperatures (50 °C and more.)



storable up to three years (50 °C and more)

What is SILPEROX

The active ingredient hydrogen peroxide breaks down into water and oxygen.




Thus,


no toxic components are available following the reaction,
fully biodegradable,
non-residual.

So, it is the most ecological and environment friend disinfectant

What is SILPEROX

CHEMISTRY





The difference from our competitors:

The others cannot maintain the level of concentration of subcomponents (hydrogen peroxide and colloidal silver) as much as what we can.

So, SILPEROX is much more stable ,more effective than the other hydrogen peroxide & colloidal silver based products.

TEST RESULTS						
ACTIVITY EVALUATION METHOD	Serial dilutions of reference Poliovirus type 1, chat strain were inoculated onto HEP-2 cells and viral titer was calculated using Spearman-Kärper method based on the virus dilution which exerted visible cytopathic effect under invert microscope.					
RESULTS		Reference virus	Effect of Silperox			
	Virus titer*	5.5	5 minute		60 minute	
			Clean Condition	Dirty Condition	Clean Condition	Dirty Condition
	Virus titer with the disinfectant **		1.5	1.5	1.5	1.5
	Reduction rate in virus titer ***		4.0	4.0	4.0	4.0
	<p>* Logarithmic TCID50 value of virus per ml</p> <p>** Logarithmic TCID50 value of virus which was exposed to disinfectant at different contact time and conditions</p> <p>*** Logarithmic TCID50 ratio of virus titer and virus titer with the disinfectant</p>					
COMMENT	<p>As tested concentrations of Silperox suspension, 10% and 1%, were found to exert cytotoxicity against test cell culture, the highest concentration of the disinfectant which did not display any toxicity, 0.1%, was used in the experiments. According to the calculations based on test results, 1 % concentration of Silperox disinfectant provided at least 4 log reduction in virus titer in all test conditions (see result table) for 5 and 60 minutes contact time. According to the Antimicrobial Division US EPA standards, disinfectants has to provide minimum 4 log virus titer reduction to be an acceptable virucidal agent.</p> <p>The results of the test show that Silperox disinfectant possesses 99.99% antiviral activity against Poliovirus type 1 at for 5 and 60 minutes contact time when used at 1% concentration.</p> <p>It has been accepted that the product, investigated for its antiviral activity against Poliovirus Type 1 which is a relevant RNA model virus according to the TS EN 14476 (March 2007) standard of Turkish Standards Institute (TSE), can display same virucidal effect against practically non-testable enveloped and non-enveloped RNA viruses such as HCV, HIV in case of using one of washing, wiping, impregnation (wetting/immersion), spraying methods at least concentration and contact time defined above.</p>					

Ayla Burçin ASUTAY
Analyst

Prof. Fikrettin ŞAHİN
Chair of Biocidal Laboratory



SILPEROX on 2019-nCoV



% 1 diluted SILPEROX is effective
% 99,99 on all kind of RNA based viruses

